

## SURGICAL CAMERA, CAMERA COUPLER, AND CAMERA CONTROL UNIT REPAIRS

### WHAT WE DO

- ◆ Cable assembly replacements.
- ◆ New glass seals.
- ◆ Prism / CCD rebuilds.
- ◆ Video boards repaired to component level.
- ◆ Camera coupler rebuilds.
- ◆ Camera control unit repairs.
- ◆ Complete testing to manufacturer specifications.

### WE REPAIR

- ◆ Surgical camera heads.
- ◆ Camera couplers.
- ◆ Surgical camera control units.

At **Jensen Instruments Technologies**, we offer surgical camera, camera coupler and camera control unit repairs for all major brands.

Our technicians are professionally trained and have many years of experience executing all major and minor repairs bringing them to like new condition.

Our repair service will only service what is necessary to return your surgical camera, camera coupler or camera control unit to its camera original condition therefore making it a cost saving venture.

We perform all major and minor repairs on surgical cameras, camera couplers and surgical camera control units manufactured by: **Olympus, Stryker, Smith & Nephew / Dyonics, Storz, Linvatec**, and many more.



Our exclusive surgical instrument repair service, Resmart Repair, promises long-term financial savings. When choosing our sustainable repair option, we offer complete product replacement if we are unable to repair the instrument to the original state.

**Repair, Refurbish, Replace.**